

Code No.16-2162

SPECIFICATION

Polyoxyethylene(20) Sorbitan Monooleate for Molecular Biology

REQUIREMENT	SPECIFICATION
Appearance	Yellow ~ yellowish brown, clear ~ slightly muddy liquid
Viscosity (25°C)	370~510mPa·s
Solubility in water	to pass test
Water	max.4.0%
Acid value	max.3
Saponification value	40~60
Hydroxyl value	65~80
Deoxyribonuclease	to pass test
Ribonuclease	to pass test