

Code No.16-2267

# SPECIFICATION

## 1-Propanol

for Molecular Biology

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Solubility in water	to pass test
Density (20°C)	0.801~0.806g/ml
Refractive index 20°C	1.383~1.387
Water	max.0.2%
Acidity (as C <sub>2</sub> H <sub>5</sub> COOH)	max.0.002%
Substances darkened by sulfuric acid	to pass test
Deoxyribonuclease	to pass test
Ribonuclease	to pass test
Assay (capillary GC)	min.99.5%