

Code No.16-2494

SPECIFICATION

Polyethylene Glycol 4,000,000

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystalline powder ~
Viscosity (0.5%, 25°C)	powder or granules
Solubility in water	under investigation
Loss on drying at 105°C	to pass test
Residue after ignition (as sulfate)	max.2.0%
pH of a 5g/l solution at 25°C	max.2.0%
Acidity (as CH ₃ COOH)	6.0~8.0
	max.0.2%