

SPECIFICATION

Potassium Dichromate

TraceSure(R)

REQUIREMENT	SPECIFICATION
* Appearance Assay (K ₂ Cr ₂ O ₇) Water-insoluble matter Chloride (Cl) Sulfate (SO ₄) Sodium (Na) Copper (Cu) Calcium (Ca) Lead (Pb) Iron (Fe)	Reddish yellow ~ yellowish red, crystals ~ crystalline powder min.99.98%(mass/mass) max.0.005%(mass/mass) max.5mg/kg max.0.03g/kg max.0.1g/kg max.1mg/kg max.0.02g/kg max.1mg/kg max.2mg/kg

JIS K 8005:2014*:Additional test performed by Wako but not required by JIS