

## SPECIFICATION

### 6 Phenols Mixture Standard Solution (Acetone Solution)

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Impurities	to pass test
2-Chlorophenol (20°C)	950~1050mg/L
Phenol (20°C)	950~1050mg/L
2,6-Dichlorophenol (20°C)	950~1050mg/L
2,4-Dichlorophenol (20°C)	950~1050mg/L
2,4,6-Trichlorophenol (20°C)	950~1050mg/L
4-Chlorophenol (20°C)	950~1050mg/L