

Code No.16-2978

SPECIFICATION

Potassium Cyanide

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder or mass
Solubility in water	to pass test
Chloride (Cl)	max.0.5%
Sulfate (SO ₄)	max.0.05%
Sulfide (S)	max.0.01%
Thiocyanate (SCN)	max.0.02%
Sodium (Na)	max.1.0%
Lead (Pb)	max.0.005%
Iron (Fe)	max.0.03%
Assay	min.95.0%