

Code No.18-0029

## SPECIFICATION

p-Rosolic Acid  
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Reddish brown, powder or mass
Solubility in ethanol	to pass test
Loss on drying at 105°C	max.15.0%(mass/mass)
Residue after ignition (as sulfate)	max.0.3%(mass/mass)
Visual transition interval (pH)	to pass test
Sensitivity	to pass test

JIS K 9037 : 2020(\*:Additional test performed by Wako but not required by JIS)