

## SPECIFICATION

## Rosmarinic Acid

for Cellbiology

REQUIREMENT	SPECIFICATION	
Appearance	White ~ pale yellowish brown, crystalline powder ~ powder or mass to pass test	
Solubility in ethanol		
Specific rotation [ $\alpha$ ]D <sub>20</sub> (c=0.2,C <sub>2</sub> H <sub>5</sub> OH) (calculated on the dehydrous basis)		+100~+110°
Water		max.3.0%
Assay (HPLC)	min.96.0%	