

Code No.18-0341

SPECIFICATION

Rutin Standard  
for Food Additives Test

REQUIREMENT	SPECIFICATION
Description	Pale yellow ~ pale yellowish green,crystalline powder to pass test min.290 - to pass test to pass test 6.9~9.1%
Identification	
Specific absorbance E1%1cm (λ max about 350nm)	
Purity	
Related substances	
* Solubility in methanol	
* Water	

This product is suited for Rutin for assay of reagents and test solutions on Japan's Specifications and Standards for Food Additives(\*:Additional test performed by Wako)