

Code No.19-0035

SPECIFICATION

Selenium Dioxide

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
Nitrate	to pass test (NO ₃ :max.about 0.01%)
Selenate and sulfate	to pass test
Copper (Cu)	max.0.002%
Lead (Pb)	max.0.003%
Tellurium (Te)	max.0.03%
Iron (Fe)	max.0.005%
Assay	min.97.0%