

Code No.19-0092

## SPECIFICATION

## Silver(I) Sulfate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White crystalline powder
Assay (Ag <sub>2</sub> SO <sub>4</sub> )	min.99.5%(mass/mass)
Solubility in dil.nitric acid	to pass test
Nitrogen compounds (as N)	max.0.001%(mass/mass)
Substances not precipitated by HCl	max.0.03%(mass/mass)
Copper (Cu)	max.0.004%(mass/mass)
Iron (Fe)	max.0.004%(mass/mass)

JIS K 8965 : 2019(\*:Additional test performed by Wako but not required by JIS)