

## SPECIFICATION

## Sodium Cyanide

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (NaCN) Solubility in water Chloride (Cl) Phosphate (PO <sub>4</sub> ) Sulfate (SO <sub>4</sub> ) Sulfide (S) Thiocyanate (SCN) Potassium (K) Lead (Pb) Iron (Fe)	White ~ slightly yellow, crystals ~ powder or plate min.97.0%(mass/mass)  to pass test  max.0.02%(mass/mass)  max.0.02%(mass/mass)  max.0.01%(mass/mass)  max.0.003%(mass/mass)  max.0.02%(mass/mass)  max.0.1%(mass/mass)  max.0.002%(mass/mass)  max.0.003%(mass/mass)

JIS K 8447 : 2018

\*:Additional test performed by Wako but not required by JIS