

SPECIFICATION

Sodium Cyanide

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (NaCN) Solubility in water Chloride (Cl) Phosphate (PO ₄) Sulfate (SO ₄) Sulfide (S) Thiocyanate (SCN) Potassium (K) Lead (Pb) Iron (Fe)	White ~ slightly yellow, crystals ~ powder or plate min.97.0%(mass/mass) to pass test max.0.02%(mass/mass) max.0.02%(mass/mass) max.0.01%(mass/mass) max.0.003%(mass/mass) max.0.02%(mass/mass) max.0.1%(mass/mass) max.0.002%(mass/mass) max.0.003%(mass/mass)

JIS K 8447 : 2018

*:Additional test performed by Wako but not required by JIS