

SPECIFICATION

Sodium Cyanide

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystals ~ powder or plate to pass test
Solubility in water	
Chloride (Cl)	max.0.1%
Sulfate (SO ₄)	max.0.1%
Sulfide (S)	max.0.02%
Thiocyanate (SCN)	max.0.02%
Lead (Pb)	max.0.005%
Iron (Fe)	max.0.03%
Assay	min.95.0%