

Code No.19-0222

## SPECIFICATION

## Sodium Phosphinate Monohydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ powder
Solubility in water	to pass test
Chloride (Cl)	max.0.01%
Sulfate (SO <sub>4</sub> )	max.0.01%
Calcium (Ca)	max.0.01%
Lead (Pb)	max.5ppm
Arsenic (As)	max.1ppm
Iron (Fe)	max.5ppm
Assay (as NaPH <sub>2</sub> O <sub>2</sub> )	82.0~86.5%

W1-4629/1997.6.5