

## SPECIFICATION

## Disodium Hydrogenphosphate Dodecahydrate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (Na <sub>2</sub> HPO <sub>4</sub> · 12H <sub>2</sub> O) Solubility in water pH of a 50g/L solution at 25°C Chloride (Cl) Nitrate Sulfate (SO <sub>4</sub> ) Potassium (K) Lead (Pb) Arsenic (As) Iron (Fe) Ammonium (NH <sub>4</sub> )	White crystals min.99.0%(mass/mass) to pass test 9.0~9.4 max.0.001%(mass/mass) to pass test max.0.005%(mass/mass) max.0.01%(mass/mass) max.5ppm(mass/mass) max.1ppm(mass/mass) max.5ppm(mass/mass) max.0.001%(mass/mass)

JIS K 9019:2021(\*:Additional test performed by Wako but not required by JIS)