

SPECIFICATION

Disodium Hydrogenphosphate Dodecahydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay (Na ₂ HPO ₄ · 12H ₂ O)	min.99.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	9.0~9.4
Chloride (Cl)	max.0.001%(mass/mass)
Nitrate	to pass test
Sulfate (SO ₄)	max.0.005%(mass/mass)
Heavy metals (as Pb)	max.5ppm(mass/mass)
Potassium (K)	max.0.01%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.5ppm(mass/mass)
Ammonium (NH ₄)	max.0.001%(mass/mass)

JIS K 9019:2016(*:Additional test performed by Wako but not required by JIS)