

## SPECIFICATION

## Potassium Sodium (+)-Tartrate Tetrahydrate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White crystals
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	6.5~8.5
Chloride (Cl)	max.0.003%
Sulfur compounds (as SO <sub>4</sub> )	max.0.03%
Copper (Cu)	max.4ppm
Lead (Pb)	max.4ppm
Iron (Fe)	max.0.002%
Ammonium (NH <sub>4</sub> )	max.0.005%
Assay	min.99.0%