

Code No.19-0315

SPECIFICATION

Sodium Salicylate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 200g/L solution at 25 °C	6.0~8.5
Sulfite and thiosulfate	to pass test (SO ₃ :max.about 0.04%, S ₂ O ₃ :max.about 0.1%)
Chloride (Cl)	max.0.03%
Sulfate (SO ₄)	max.0.05%
Heavy metals (as Pb)	max.0.002%
Assay	min.99.0%