

## SPECIFICATION

## Sodium Sulfite

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (Na <sub>2</sub> SO <sub>3</sub> ) Solubility in water pH of a 50g/l solution at 25°C Chloride (Cl) Potassium (K) Copper (Cu) Calcium (Ca) Zinc (Zn) Lead (Pb) Arsenic (As) Iron (Fe)	White crystalline powder min.97.0%(mass/mass) to pass test 9.0~11.0 max.0.02%(mass/mass) max.0.05%(mass/mass) max.5ppm(mass/mass) max.0.01%(mass/mass) max.5ppm(mass/mass) max.0.001%(mass/mass) max.1ppm(mass/mass) max.0.001%(mass/mass)

JIS K 8061:2010(\*:Additional test performed by Wako but not required by JIS)