

SPECIFICATION

Sulfanilic Acid

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|---|---|
| * Appearance Assay(C ₆ H ₇ NO ₃ S) Solubility in sodium carbonate solution Residue after ignition (as sulfate) Chloride (Cl) Sulfate (SO ₄) Nitrite (NO ₂) | White ~ slightly brown, crystalline powder min.99.0%(mass/mass) to pass test max.0.02%(mass/mass) max.0.002%(mass/mass) max.0.005%(mass/mass) max.0.5ppm(mass/mass) |

JIS K 8586:2011(*:Additional test performed by Wako but not required by JIS)