

Code No.19-0471

## SPECIFICATION

## Sulfuric Acid, Fuming

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	to pass test
Assay (free SO <sub>3</sub> )	60.0~65.0%(mass/mass)
Residue after ignition (as sulfate)	max.0.002%(mass/mass)
Chloride (Cl)	max.2ppm(mass/mass)
Nitrate	to pass test
Heavy metals (as Pb)	max.3ppm(mass/mass)
Arsenic (As)	max.0.1 ppm(mass/mass)
Selenium (Se)	max.5ppm(mass/mass)
Iron (Fe)	max.4ppm(mass/mass)
Substances reducing permanganate	to pass test

JIS K 8741:2018

\*:Additional test performed by Wako but not required by JIS  
The assay of this product is approximately 60%