

Code No.19-0967

SPECIFICATION

6mol/l Sodium Thiocyanate Solution for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Chloride (Cl)	max.0.005%
Sulfate (SO ₄)	max.0.003%
Copper (Cu)	max.1ppm
Mercury (Hg)	max.0.5ppm
Lead (Pb)	max.1ppm
Iron (Fe)	max.1ppm
Ammonium (NH ₄)	max.0.003%
Iodine consuming substances (as S)	max.0.004%
Filtration test	to pass test
Concentration (20°C)	5.8~6.2mol/l