

Code No.19-1003

# SPECIFICATION

## Sodium Taurocholate

Wako Special Grade

| REQUIREMENT                                | SPECIFICATION  |
|--|--|
| Appearance                                 | White ~ slightly brown,<br>crystalline powder ~ powder or mass<br>to pass test |
| Solubility in water                        |  |
| Specific rotation [ $\alpha$ ]D20(c=3,H2O) | +20~+24°   |
| Water                                      | max.10.0%  |
| pH of a 30g/l solution at 25 °C            | 4.5~7.5  |
| Thin layer chromatography test             | to pass test   |