

Code No.19-1952

SPECIFICATION

Silicone Resin

for the Japanese Pharmacopoeia General Tests

REQUIREMENT	SPECIFICATION
Description	Colorless ~ light gray, viscous liquid or a pasty material
Refractive index n _{D25}	1.400~1.410
Kinematic viscosity (25°C)	100~1100mm ² /s
Hexane insoluble substances	not more than 2.25g
Specific gravity d ₂₀₂₀	0.96~1.02

This product is suited for Silicone Resin of general tests on The Japanese Pharmacopoeia.