

SPECIFICATION

Mercaptoacetic Acid

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (HSCH ₂ COOH) Solubility in water Residue after ignition (as sulfate) Suitability for iron determination Iron (Fe)	Colorless ~ nearly colorless, clear liquid min.90.0%(mass/mass) to pass test max.0.01%(mass/mass) to pass test max.5ppm(mass/mass)

JIS K 8630:2019(*:Additional test performed by Wako but not required by JIS)