

Code No.20-0145

## SPECIFICATION

## Thymolphthalein

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly yellow, crystalline powder ~ powder to pass test
Solubility in ethanol	
Absorbance (5mg/l, pH10.5)	min.0.35
Loss on drying at 105°C	max.1.0%(mass/mass)
Residue after ignition (as sulfate)	max.0.1%(mass/mass)
Visual transition interval (pH)	8.6~10.5 (colorless ~ blue)

JIS K 8642:2011

(\*:Additional test performed by Wako but not required by JIS)