

Code No.20-0157

SPECIFICATION

Tin(II) Chloride Dihydrate

Wako 1st Grade

| REQUIREMENT | SPECIFICATION |
|-------------------------------------|----------------|
| Appearance | White crystals |
| Solubility in dil.hydrochloric acid | to pass test |
| Sulfate (SO ₄) | max.0.01% |
| Copper (Cu) | max.0.002% |
| Lead (Pb) | max.0.005% |
| Iron (Fe) | max.0.005% |
| Assay | min.95.0% |