

SPECIFICATION

Toluene

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (C ₆ H ₅ CH ₃) (GC) Density (20°C) Refractive index 20°C Water Residue after evaporation Acidity (as HCl) Alkalinity (as NaOH) Sulfur compounds (as S) Thiophene Substances darkened by sulfuric acid	Colorless clear liquid min.99.5%(mass/mass) 0.864~0.868g/ml 1.496~1.498 max.0.03%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.003%(mass/mass) to pass test to pass test

JIS K 8680:2006*:Additional test performed by Wako but not required by JIS