

Code No.20-0348

SPECIFICATION

Turpentine Oil

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	0.853~0.863g/ml
Refractive index 20°C	1.463~1.470
Specific rotation [α]D20(c=10,C2H5OH)	-15~-35°
Substances darkened by hydrochloric acid	to pass test
Substances darkened by potassium hydroxide	to pass test
α -Pinene + β -Pinene (capillary GC)	min.90.0%