

SPECIFICATION

Turpentine Oil

Wako 1st Grade

| REQUIREMENT | SPECIFICATION |
|---|---|
| Appearance | Colorless ~ slightly yellow, clear liquid |
| Density (20°C) | 0.853~0.863g/ml |
| Refractive index 20°C | 1.463~1.470 |
| Specific rotation [α]D20(c=10,C2H5OH) | -15~-35° |
| Substances darkened by hydrochloric acid | to pass test |
| Substances darkened by potassium hydroxide | to pass test |
| α -Pinene + β -Pinene (capillary GC) | min.90.0% |