

SPECIFICATION

2-Methyltetrahydrofuran

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ yellowish brown, clear liquid
Density (20°C)	0.850~0.860g/ml
Refractive index 20°C	1.402~1.409
Water	max.0.1%
Peroxides (as H ₂ O ₂)	max.0.2%
Assay (capillary GC)	min.98.0%
2,6-Di-t-butyl-4-methylphenol (stabilizer)	0.01~0.05%