

Code No.20-0631

SPECIFICATION

N,N,N',N'-Tetramethylethylenediamine  
for Electrophoresis

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Solubility in water	to pass test
Density (20°C)	0.772~0.778g/mL
Refractive index 20°C	under investigation
Water	under investigation
Residue after evaporation	max.0.05%
Heavy metals (as Pb)	max.0.001%
Iron (Fe)	max.0.001%
Assay (capillary GC)	min.99.0%