

Code No.20-0706

SPECIFICATION

1,1,2,2-Tetrachloroethane

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|--|
| * Appearance | Colorless ~ nearly colorless, clear liquid |
| Assay (CHCl ₂ CHCl ₂)(GC) | min.97.0%(mass/mass) |
| Density (20°C) | 1.593~1.601g/ml |
| Refractive index 20°C | 1.493~1.497 |
| Water | max.0.1%(mass/mass) |
| Acidity (as HCl) | max.0.003%(mass/mass) |

JIS K 9520:2013

(*:Additional test performed by Wako but not required by JIS)