

SPECIFICATION

1 mol/L Triethylammonium Hydrogencarbonate Solution for Nucleic Acid Synthesis

| REQUIREMENT | SPECIFICATION |
|--------------------|------------------------|
| Appearance | Colorless clear liquid |
| pH at 15°C | 7.3~7.8 |
| Absorbance (220nm) | 0.3~0.9 |
| Absorbance (230nm) | 0.1~0.5 |
| Absorbance (270nm) | max.0.05 |
| Concentration | 0.95~1.05mol/L |