

Code No.20-0840

SPECIFICATION

 $(\pm)$ - α -Tocopherol Acetate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ yellow, clear liquid
Refractive index 20°C	1.494~1.499
Specific absorbance E1%1cm (284nm)	41.0~45.0
Water	max.0.2%
Acidity (as CH ₃ COOH)	max.0.2%
α -Tocopherol	max.0.5%
Assay (HPLC)	min.98.0%