

## SPECIFICATION

## Tetrahydrofuran

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (C <sub>4</sub> H <sub>8</sub> O) (GC) Density (20°C) Refractive index 20°C Water Peroxides (as H <sub>2</sub> O <sub>2</sub> )	Colorless clear liquid min.99.5%(mass/mass) 0.884~0.889g/ml 1.406~1.409 max.0.05%(mass/mass) max.0.005%(mass/mass)

JIS K 9705:2013\*:Additional test performed by Wako but not required by JIS(with antioxidants)