

SPECIFICATION

2-Ethyl-2-hydroxymethyl-1,3-propanediol Trimethacrylate

Practical Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellow, clear liquid
Solubility in ethanol	under investigation
Density (20°C)	1.062~1.072g/ml
Refractive index 20°C	1.470~1.476
Water	max.0.2%
Acidity [as CH ₂ :C(CH ₃)COOH]	max.0.2%
Assay (capillary GC)	min.90.0%