

SPECIFICATION

Titanium(III) Chloride Solution

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (TiCl ₃) Nitrogen compounds (as N) Magnesium (Mg) Lead (Pb) Iron (Fe) Nickel (Ni)	Dark purple liquid min.20%(mass/mass) max.0.01%(mass/mass) max.0.005%(mass/mass) max.0.002%(mass/mass) max.0.005%(mass/mass) max.0.002%(mass/mass)

JIS K 8401:2011

(*:Additional test performed by Wako but not required by JIS)