

## SPECIFICATION

## Tetrahydrofuran, Stabilizer Free

## Infinity Pure

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.886~0.890g/mL
Refractive index 20°C	1.406~1.409
Absorbance (230nm)	max.0.50
Water	max.0.02%
Residue after evaporation	max.10ppm
Acidity (as CH <sub>3</sub> COOH)	max.10ppm
Copper (Cu)	max.5ppb
Lead (Pb)	max.5ppb
Iron (Fe)	max.10ppb
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.0.01%
Fluorescence test	to pass test
Assay (capillary GC)	min.99.8%