

SPECIFICATION

Tetrahydrofuran, Stabilizer Free

Infinity Pure

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.886~0.890g/ml
Refractive index 20°C	1.406~1.409
Absorbance (230nm)	max.0.50
Water	max.0.02%
Residue after evaporation	max.10ppm
Acidity (as CH ₃ COOH)	max.10ppm
Copper (Cu)	max.5ppb
Lead (Pb)	max.5ppb
Iron (Fe)	max.10ppb
Peroxides (as H ₂ O ₂)	max.0.01%
Fluorescence test	to pass test
Assay (capillary GC)	min.99.8%