

Code No.20-1604

SPECIFICATION

Trometamol

Japanese Pharmaceutical Codex

REQUIREMENT	SPECIFICATION
* Description	White crystalline powder
Identification	to pass test
pH	10.3~10.7
Melting point	168~172°C
Purity	-
(1) Clarity and color of solution	to pass test
(2) Heavy metals	max.8ppm
(3) Arsenic	max.1.6ppm
Loss on drying	max.0.5%
Residue on ignition	max.0.10%
Assay (after drying)	min.99.0%

Japanese Pharmaceutical Codex 2002(*:Reference item)