

Code No.20-1700

SPECIFICATION

Tetramethylammonium Chloride

for Molecular Biology

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystalline powder
Solubility in water	to pass test
Solubility in ethanol	to pass test
Water	max.2.0%
pH of a 50g/l solution at 25°C	4.5~6.5
Sulfate (SO ₄)	max.0.005%
Iron (Fe)	max.0.001%
Deoxyribonuclease	to pass test
Ribonuclease	to pass test
Assay	min.98.0%