

## SPECIFICATION

Tetrahydrofuran, Super Dehydrated, Stabilizer Free  
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.884~0.889g/mL
Refractive index 20°C	1.406~1.409
Absorbance (250nm)	max.0.50
Absorbance (280nm)	max.0.20
Absorbance (350nm)	max.0.02
Water	max.0.003%
Acidity (as CH <sub>3</sub> COOH)	max.0.002%
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.0.02%
Assay (capillary GC)	min.99.5%