

Code No.20-1856

## SPECIFICATION

## Trichloroethylene

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (CICH=CCl <sub>2</sub> )(GC)	min.99.5%(mass/mass)
Density (20°C)	1.461~1.469g/ml
Refractive index 20°C	1.477~1.479
Water	max.0.01%(mass/mass)
Residue after evaporation	max.0.001%(mass/mass)
Acidity (as HCl)	max.0.001%(mass/mass)
Substances reducing ammoniac silver nitrate	to pass test
Free chlorine	to pass test
Substances darkened by sulfuric acid	to pass test

JIS K 8666:2013(\*:Additional test performed by Wako but not required by JIS)