

## SPECIFICATION

### Transforming Growth Factor- $\beta$ 3, Human, recombinant, Animal-derived-free for Cellbiology

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Bacterial endotoxins	less than 0.01 ng/ $\mu$ g
Electrophoresis test (SDS-PAGE)	to pass test
Mass spectrometry (MALDI-TOF)	to pass test

Source: E. coli expressed human transforming growth factor- $\beta$  3 This product is recombinant cytokine which is produced in the process of culture and purification without using any substances from animals.