

Code No.20-1977

SPECIFICATION

Transferrin(Apo), from Human  
for Biochemistry

| REQUIREMENT                     | SPECIFICATION                         |
|---------------------------------|---------------------------------------|
| Appearance                      | White ~ pale brown, crystals ~ powder |
| Solubility in water             | to pass test                          |
| Loss on drying at 105°C         | under investigation                   |
| Iron (Fe)                       | max.0.05%                             |
| Electrophoresis test (SDS-PAGE) | to pass test                          |
| Assay (elemental analysis)      | 98~105%                               |