

SPECIFICATION

Tetrahydrofuran, with Stabilizer
for GPC

| REQUIREMENT | SPECIFICATION |
|---|------------------------|
| Appearance | Colorless clear liquid |
| Density (20°C) | 0.884~0.889g/ml |
| Refractive index 20°C | 1.406~1.409 |
| Absorbance (300nm) | max.1.00 |
| Absorbance (320~600nm) | max.0.10 |
| Water | max.0.05% |
| Peroxides (as H ₂ O ₂) | max.0.005% |
| Suitability for GPC test | to pass test |
| Assay (capillary GC) | min.99.5% |
| Stabilizer (BHT) | 0.01~0.03% |