

SPECIFICATION

Uranine

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Solubility in water Ethanol-insoluble matter Absorbance (3mg/l, pH7.0) Loss on drying at 105°C Residue after ignition (as sulfate)(after drying) Chloride (Cl) Zinc (Zn) Sensitivity	Green or yellowish red ~ reddish brown, crystalline powder ~ powder or mass to pass test max.1.0% min.0.45 max.7.0% 35.0~40.5% max.0.1% max.0.05% to pass test

JIS K 8830-1995

(*:Additional test performed by Wako but not required by JIS)