

SPECIFICATION

Ultrapure Water
for Ultratrace Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Lithium (Li)	max.10ppt
Sodium (Na)	max.10ppt
Potassium (K)	max.10ppt
Rubidium (Rb)	max.10ppt
Copper (Cu)	max.10ppt
Silver (Ag)	max.10ppt
Beryllium (Be)	max.10ppt
Magnesium (Mg)	max.10ppt
Calcium (Ca)	max.10ppt
Strontium (Sr)	max.10ppt
Barium (Ba)	max.10ppt
Zinc (Zn)	max.10ppt
Cadmium (Cd)	max.10ppt
Aluminium (Al)	max.10ppt
Thorium (Th)	max.10ppt
Tin (Sn)	max.10ppt
Lead (Pb)	max.10ppt
Arsenic (As)	max.10ppt
Bismuth (Bi)	max.10ppt
Chromium (Cr)	max.10ppt
Molybdenum (Mo)	max.10ppt
Tungsten (W)	max.10ppt
Manganese (Mn)	max.10ppt
Iron (Fe)	max.10ppt
Cobalt (Co)	max.10ppt
Rhodium (Rh)	max.10ppt
Nickel (Ni)	max.10ppt
Uranium (U)	max.10ppt

